Introduction to sequencing: Sequencing technology and preprocessing of sequencing data

Koen Van den Berge

Last compiled on 09 November, 2021

Contents

1 The Smart-seq2 protocol

1

1 The Smart-seq2 protocol

- 1. First, the biological **sample of interest are collected**, and cells should be disaggregated into a single-cell suspension.
- 2. Cells are picked from the suspension and embedded in a **lysis** buffer containing a ribonuclease inhibitor that blocks RNA degradation, and primers for the reverse transcription step.
- 3. Reverse transcription. Single-stranded RNA molecules are converted into more stable double-stranded cDNA molecules.
- 4. The cDNA molecules are amplified using **PCR amplification**.
- 5. Construction of sequencing libraries using **tagmentation**. A transposase is used to **fragment** the cDNA molecules and simultaneously **ligate adapter sequences**.
- 6. Another cycle of PCR occurs (the enrichment PCR step). Sample barcodes may be added to allow multiplexing; a a dual-index strategy is often used, indices referred to as index 1 (i7) and index 2 (i5) in the Figure below.

1.0.1 The 10X Genomics protocol




Figure 1: Figure: The SMART-Seq2 library prep workflow. Image from Picelli $et\ al.\ (2014).$